



Be Ready. When Minutes Matter[®]

Pathogen Reduced Cryoprecipitated Fibrinogen Complex
(INTERCEPT[®] Fibrinogen Complex)

produced from the INTERCEPT[®] Blood System for Cryoprecipitation

INTERCEPT Fibrinogen Complex is available for immediate use for up to 5 days when stored thawed; and when stored frozen requires thawing prior to use.

INTERCEPT® Fibrinogen Complex


Breakthrough Device to CONTROL BLEEDING

The ready-to-use INTERCEPT® Fibrinogen Complex (IFC)* is approved specifically for the treatment and control of bleeding, including massive hemorrhage, associated with fibrinogen deficiency.

- Immediate, enriched source of key factors in effective hemostasis¹⁻³
 - » Fibrinogen
 - » Factor XIII
 - » von Willebrand Factor
 - » Other vital clotting proteins




5 Day Room Temperature Shelf Life



TRANSFUSE IMMEDIATELY With the First Blood Components*

Thaw IFC in advance to minimize wait times and have available for immediate use.



* INTERCEPT Fibrinogen Complex is available for immediate use for up to 5 days when stored thawed; and when stored frozen requires thawing prior to use.
**There is no pathogen inactivation process that has been shown to eliminate all pathogens. Certain non-enveloped viruses (e.g., hepatitis A virus (HAV), hepatitis E virus (HEV), parvovirus B19 and poliovirus) and Bacillus cereus spores have demonstrated resistance to the INTERCEPT process.


PROTECT PATIENTS: Pathogen Reduction

IFC provides broad spectrum transfusion transmitted infection (TTI) risk reduction, including viruses, bacteria, and emerging pathogens.^{8,9,**}

- Produced from plasma treated by the INTERCEPT Blood System.
- Pathogen reduction enables IFC’s 5-day post-thaw shelf life.

The INTERCEPT® Blood System has 20 years of clinical and post-market surveillance experience.

INTERCEPT® Blood System for Plasma Mechanism of Action




Upon UVA illumination, amotosalen cross-links nucleic acids to block replication and inactivates pathogens

IMPROVE EFFICIENCIES: Turnaround Times and Wastage


IFC accelerates availability, minimizes wait times and increases reliability and predictability of blood component delivery.

Improved Turnaround Times



IFC reduces blood component wastage, improving blood stewardship.

Significant Wastage Reduction



The INTERCEPT® Fibrinogen Complex Advantage

TRANSFUSE
IMMEDIATELY*



IMPROVE
EFFICIENCIES



CONTROL
BLEEDING†



PROTECT
PATIENTS‡



*INTERCEPT Fibrinogen Complex is available for immediate use for up to 5 days when stored thawed; and when stored frozen requires thawing prior to use.

†Bleeding associated with fibrinogen deficiency.

‡Broad spectrum transfusion transmitted infection risk reduction.

Availability

Ready for Immediate Use!

Once thawed, may be stored at room temperature for up to 5 days.

- Provided in single-use containers
- Components may be purchased as single or pre-pooled units



Catalog #	Description	Average Fibrinogen (mg)*	# Donors†
FC10	Pooled Fibrinogen Complex 1.0, Cryoprecipitated, Psoralen Treated	740	2
FC15	Pooled Fibrinogen Complex 1.5, Cryoprecipitated, Psoralen Treated	1,556	4
FC20	Pooled Fibrinogen Complex 2.0, Cryoprecipitated, Psoralen Treated	2,220**	6
FC30	Pooled Fibrinogen Complex 3.0, Cryoprecipitated, Psoralen Treated	2,845	8
FC40	Pooled Fibrinogen Complex 4.0, Cryoprecipitated, Psoralen Treated	3,700**	10

* Mean fibrinogen content (per package insert). Fibrinogen content depends on donor plasma fibrinogen levels.

** Calculated based on pooling of FC10.

† Number of donors based on whole blood donors.

Hemorrhage Increases Mortality In



Hemorrhage is a Leading Cause of Preventable Death²³

Trauma is the

#1

cause of death

in adults <45 years old²⁴

Of trauma deaths

~40%

from hemorrhage²⁴

Hours until death

~1.6

from exsanguination²⁵

Faster is Better in Hemorrhage Control²⁶

Massive Transfusion Protocols (MTP) were developed to improve hemorrhage outcomes by delivering blood components quickly.

- Every minute of delay between the activation of an MTP and the arrival of the first blood components, results in a 5% increase in the odds of mortality.
- Timely delivery of blood components is an important metric, similar to “door-to-balloon” time.

EVERY
1
MINUTE

=

5%
RISK OF
DEATH


DELAY INCREASES MORTALITY²⁶

*Severe loss of blood

Effective Treatment: Restoring Fibrinogen & Other Clotting Factors


Early fibrinogen supplementation restores clot strength, reduces blood loss, and lowers mortality¹⁸

- Fibrinogen and other key clotting factors decrease rapidly and significantly in hemorrhage.^{18,27}
- Fibrinogen level is an independent risk factor for hemorrhage and mortality in trauma,²⁸ cardiac (CV) surgery²⁹ and postpartum hemorrhage.³⁰



Early delivery of **Fibrinogen**, **factor XIII** and **von Willebrand factor** adds the clotting strength needed to achieve stable clot formation and restore hemostasis.^{18,31,32}

Early Fibrinogen Supplementation Improves Patient Outcomes*



*These results are based on studies evaluating cryo AHF. Results with IFC may vary.

MTPs Lack Critical Components from the Start

Cryoprecipitated AHF Inventory Challenges

LONG WAIT TIMES
15 min – 2.8 hours^{4,7}
Stored frozen and typically thawed in round 2 or 3 of MTP³⁶


SHORT SHELF LIFE
4 – 6 hours³⁸
Post thaw due in part to infectious risk

HIGH WASTAGE RATES
7 – 33%³⁷
Thawed cryo AHF is wasted

In >75% of U.S. exsanguination cases, cryo AHF arrives too late to be medically efficacious³⁶

Delays Impact Cryoprecipitated AHF Utility

Time to issue cryoprecipitated AHF by patient type^{39*}



- **2.5 hours:** Median time to receipt of cryo AHF after initiation of an MTP³⁹
 - **Trauma and OB earliest: median of 1.7 hours**³⁹
- **1.6 hours:** Median time to death from exsanguination²⁵

*©2017 John Wiley & Sons Ltd. Graph reproduced with permission.

Intended Use

• Treatment and control of bleeding, including massive hemorrhage, associated with fibrinogen deficiency. • Control of bleeding when recombinant and/or specific virally inactivated preparations of factor XIII or von Willebrand factor (vWF) are not available. • Second-line therapy for von Willebrand disease (vWD). • Control of uremic bleeding after other treatment modalities have failed.

Limitations of Use: Pathogen Reduced Cryoprecipitated Fibrinogen Complex should not be used for replacement of factor VIII.

Contraindications

Contraindicated for preparation of blood components intended for patients with a history of hypersensitivity reaction to amotosalen or other psoralens.

Contraindicated for preparation of blood components intended for neonatal patients treated with phototherapy devices that emit a peak energy wavelength less than 425 nm, or have a lower bound of the emission bandwidth <375 nm, due to the potential for erythema resulting from interaction between ultraviolet light and amotosalen.

Warnings and Precautions

Only the INTERCEPT Blood System for Cryoprecipitation is approved for use to produce Pathogen Reduced Cryoprecipitated Fibrinogen Complex.

For management of patients with vWD or factor XIII deficiency, Pathogen Reduced Cryoprecipitated Fibrinogen Complex should not be used if recombinant or specific virally-inactivated factor preparations are available. In emergent situations, if recombinant or specific virally-inactivated factor preparations are not available, Pathogen Reduced Cryoprecipitated Fibrinogen Complex may be administered.

References

1. Levy JH, Welsby I, et al. Fibrinogen as a therapeutic target for bleeding: a review of critical levels and replacement therapy. *Transfusion* 2014;54(5):1389-1405; quiz 1388.
2. Schroeder V, Kohler HP. Factor XIII: Structure and Function. *Semin Thromb Hemost* 2016;42(4):422-428.
3. Peyvandi, F. Diagnosis and management of patients with von Willebrand's disease in Italy: an Expert Meeting Report. *Blood Transfus* 2018;16(4):326-328.
4. AABB. Circular of Information for the Use of Human Blood and Blood Components. Bethesda, MD: AABB; 2024.
5. RIASTAP PI. CSL Behring LLC; 2021.
6. FIBRYGA PI. Octapharma; 2024.
7. Holcomb JB, Fox EE, Zhang X, et al. Cryoprecipitate Use in the Prospective Observational Multicenter Major Trauma Transfusion study (PROMTTT). *The journal of trauma and acute care surgery* 2013;75:S31-S39.
8. INTERCEPT Blood System for Plasma Package Insert.
9. INTERCEPT Blood System for Cryoprecipitation for the manufacturing of Pathogen Reduced Cryoprecipitated Fibrinogen Complex Package Insert.
10. Sethapati VR, Pham TD, Quach T, Nguyen A, Le J, Cai W, Virk, MS. Implementation and early outcomes with Pathogen Reduced Cryoprecipitated Fibrinogen Complex. *American Journal of Clinical Pathology* 2024;1-6.
11. Aidikoff J, Trivedi D, Kwok R, Shafi H. How do I implement pathogen reduced Cryoprecipitated fibrinogen complex in a tertiary Hospital's blood Bank. *Transfusion* 2024;1-10.
12. Webb C, Hanna D, Zhang R, Sarode R. Improved Operational Efficiencies With Pathogen Reduced Cryoprecipitated Fibrinogen Complex vs. CryoAHF. (2023), P-TS-43. *Transfusion*, 63: 283A-284A.
13. Barnes-Jewish Hospital (2024). Cerus Corporation.
14. Goertzen S, Wu S. Implementation of Pathogen Reduced Cryoprecipitated Fibrinogen Complex in a Pediatric Level II Trauma Center. (2023), P-TS-41. *Transfusion*, 63: 282A-282A.
15. Alvarado A. (2024). Cerus Corporation.
16. Cushing M, et al. (2024), P-TS-118. AABB 2024.
17. Stanworth SJ, Davenport R, Curry N, et al. Mortality from trauma haemorrhage and opportunities for improvement in transfusion practice. *The British journal of surgery* 2016;103:357-65.
18. Rourke C, Curry N, Khan S, et al. Fibrinogen levels during trauma hemorrhage, response to replacement therapy, and association with patient outcomes. *Journal of thrombosis and haemostasis:JTH* 2012;10:1342-51.
19. Butwick AJ, Goodnough LT. Transfusion and coagulation management in major obstetric hemorrhage. *Current opinion in anaesthesiology* 2015;28:275-84. Stanworth SJ, Davenport R, Curry N, et al. Mortality from trauma haemorrhage and opportunities for improvement in transfusion practice. *The British journal of surgery* 2016;103:357-65.
20. Görlinger K, Shore-Lesserson L, Dirkmann D, Hanke AA, Rahe-Meyer N, Tanaka KA. Management of hemorrhage in cardiothoracic surgery. *J Cardiothorac Vasc Anesth* 2013;27:S20-34.
21. Stinger HK, Spinella PC, Perkins JG, et al. *J Trauma*. 2008;64:S79-S85.
22. Joint Trauma System, Damage Control Resuscitation Clinical Practice Guideline, 12 July 2019. [https://jts.amedd.army.mil/assets/docs/cpgs/JTS_Clinical_Practice_Guidelines_\(CPGs\)/Damage_Control_Resuscitation_12_Jul_2019_ID18.pdf](https://jts.amedd.army.mil/assets/docs/cpgs/JTS_Clinical_Practice_Guidelines_(CPGs)/Damage_Control_Resuscitation_12_Jul_2019_ID18.pdf). Accessed Jul 2020.
23. Drake SA, Holcomb JB, Yang Y, et al. Establishing a Regional Trauma Preventable/Potentially Preventable Death Rate. *Annals of surgery* 2018.
24. Callcut RA, Kornblith LZ, Conroy AS, et al. The why and how our trauma patients die: A prospective Multicenter Western Trauma Association study. *The journal of trauma and acute care surgery* 2019;86:864-70.
25. Cripps MW, Kutcher ME, Daley A, et al. Cause and timing of death in massively transfused trauma patients. *The journal of trauma and acute care surgery* 2013;75:S255-62.
26. Meyer DE, Vincent LE, Fox EE, et al. Every minute counts: Time to delivery of initial massive transfusion cooler and its impact on mortality. *The journal of trauma and acute care surgery* 2017;83:19-24.
27. Burggraf M, Payas A, Kauther MD, Schoeneberg C, Lendemans S. Evaluation of clotting factor activities early after severe multiple trauma and their correlation with coagulation tests and clinical data. *World J Emerg Surg*. 2015 Sep 22;10:43.
28. McQuilten ZK, Wood EM, Bailey M, Cameron PA, Cooper DJ. Fibrinogen is an independent predictor of mortality in major trauma patients: A five-year statewide cohort study. *Injury* 2017;48:1074-81.
29. Ranucci M, Pistuddi V, Baryshnikova E, Colella D, Bianchi P. Fibrinogen Levels After Cardiac Surgical Procedures: Association With Postoperative Bleeding, Trigger Values, and Target Values. *The Annals of Thoracic Surgery* 2016;102:78-85.
30. Charbit B, Mandelbrot L, Samain E, et al. The decrease of fibrinogen is an early predictor of the severity of postpartum hemorrhage. *Journal of thrombosis and haemostasis: JTH* 2007;5:266-73.
31. Levy JH, Szlam F, Tanaka KA, Sniecinski RM. Fibrinogen and hemostasis: a primary hemostatic target for the management of acquired bleeding. *Anesthesia and analgesia* 2012;114:261-74.
32. Chapin JC, Hajjar KA. Fibrinolysis and the control of blood coagulation. *Blood reviews* 2015;29:17-24.
33. Dittilo M, et al. The role of cryoprecipitate in massively transfused patients: Results from the Trauma Quality Improvement Program database may change your mind. *J Trauma Acute Care Surg*. 2020. 89(2): p. 336-343.
34. Pearse BL, et al. Protocol guided bleeding management improves cardiac surgery patient outcomes. *Vox Sang*. 2015. 109(3): p. 267-79.
35. Green L, et al. Early cryoprecipitate transfusion versus standard care in severe postpartum haemorrhage: a pilot cluster-randomised trial. *Anaesthesia*. 2022. 77(2): p. 175-184.
36. Data on file. Calculation based on references: (17) Cripps et al. 2013 and (4) Holcomb et al. 2013.
37. Dunbar NM, Olson NJ, Szczepiorkowski ZM, et al. Blood component transfusion and wastage rates in the setting of massive transfusion in three regional trauma centers. *Transfusion* 2017;57:45-52.
38. Wagner SJ, Hapip CA, Abel L. Bacterial safety of extended room temperature storage of thawed crwyoprecipitate. *Transfusion* 2019;59:3549-50.
39. McQuilten ZK, Bailey M, Cameron PA, Stanworth SJ, Venardos K, Wood EM, Cooper DJ. Fibrinogen concentration and use of fibrinogen supplementation with cryoprecipitate in patients with critical bleeding receiving massive transfusion: a bi-national cohort study. *Br J Haematol*. 2017 Oct;179(1):131-141.

Rx only.



Global Headquarters | 1220 Concord Avenue | Concord, CA US 94520 | 855.835.3523

cerus.com | INTERCEPTFibrinogenComplex.com