

# Impact of Pathogen Reduction (PR) vs. LVDS Testing on Platelet Availability: A Study Based on Real-World Experience



Tracy Collier, Vera Chrebtow  
Cerius Corporation, Concord, CA

## Background

Blood shortages continue to raise interest in optimizing blood availability. Two bacterial risk mitigation strategies under the FDA are pathogen reduction (PR) and large volume delayed sampling at 48-hours (LVDS 48hr). These methods differentially impact platelet component (PC) availability for release, shelf-life, time to transfusion, and waste. An independent blood center (BC) that supplies >10,000 PC to more than 30 hospitals annually, assessed PC collection and distribution data to evaluate PC availability when comparing PR (INTERCEPT® Blood System) and LVDS 48hr. PC access also was evaluated at a level II trauma center hospital serviced by the BC. The BC and hospital requested to remain anonymous but reviewed the abstract and agreed with the results and conclusions.

## Aims

To evaluate platelet component availability when comparing pathogen reduction (INTERCEPT® Blood System) and large volume delayed sampling 48-hour.

## Methods

PC data from October 2021 through February 2022 were exported from the blood establishment computer software (BECS) including DIN, PC type, collection date/time, and shipment date/time. PC age at distribution and usable shelf-life were compared between LVDS 48hr and PR; usable shelf-life was calculated based on maximal shelf-life (7 and 5 days for LVDS 48hr and PR, respectively) and ship date/time. PC wastage at BC was also assessed. The hospital provided one month of data including the date/time of transfusion. Data was cross-referenced with the collection information from the BC; the PC age upon transfusion was determined and compared between LVDS 48hr and PR.

## Results

The BC distributed 4,793 components during the study period; 90% were PR PC while the remainder were LVDS 48hr PC (**Table 1**). Analysis of collection and distribution data demonstrated average PR PC release 64 hours (2.7 days) earlier with greater remaining usable


shelf-life compared to LVDS 48hr PC. Earlier release translated to significantly fewer wasted PC with PR (**Figure 1**). At the hospital, most of the transfusions occurred between day 3-5 (82%) vs. LVDS 48hr between day 5 and 7 (76%) (**Figure 2**).

**Table 1: Blood Center Distribution, Shelf Life, and Wastage**


	PC Type		p-value
	PR	LVDS 48Hr	
Total # Units	4326	467	--
Average Age at Distribution (Hr )	47.7	112.4	<0.0001*
Usable Shelf-Life at Distribution (Hr )	96.3	79.6	<0.0001*
# Units Wasted (%)	341 (7.9)	120 (25.7)	<0.0001**

\* T-Test; p<0.05 is statistically significant. \*\*Chi-square test; p<0.05 is statistically significant.

**Figure 1: Platelet Average Age at Release and Remaining Shelf-Life at the Blood Center**



**Figure 2: Distribution of Platelet Age at Time of Transfusion**



The hospital routinely accepts short-dated platelets thus resulting in skewed distribution/transfusion toward the end of expiry for both platelet types (i.e. Day 4-5 for PR and Day 7 for LVDS).

## Conclusions

In this study, earlier release of PR PC enabled sooner availability of PC and an extended usable shelf-life, particularly if blood center production such as infectious disease testing turn-around time is optimized. Conversely, LVDS 48hr resulted in delayed sampling and a decreased usable shelf-life with increased waste. From a hospital and patient perspective, the comparison of PC age at transfusion demonstrated earlier availability and transfusion of fresher platelets with PR.